

MAP3K11

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 93kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human MAP3K11

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MLK3; PTK1; SPRK; MLK-3; MEKK11;

Catalog #: A7623

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4296

Isotype: IgG

Swiss Prot: Q16584

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the serine/threonine kinase family. This kinase contains a SH3 domain and a leucine zipper-basic motif. This kinase preferentially activates MAPK8/JNK kinase, and functions as a positive regulator of JNK signaling pathway. This kinase can directly phosphorylate, and activates I κ B kinase alpha and beta, and is found to be involved in the transcription activity of NF- κ B mediated by Rho family GTPases and CDC42.

To place an order, please [Click HERE](#).