

PKP3

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 87kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human PKP3

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A7624

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 11187

Isotype: IgG

Swiss Prot: Q9Y446

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the arm-repeat (armadillo) and plakophilin gene families. Plakophilin proteins contain numerous armadillo repeats, localize to cell desmosomes and nuclei, and participate in linking cadherins to intermediate filaments in the cytoskeleton. This protein may act in cellular desmosome-dependent adhesion and signaling pathways. Two transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).