

CAST

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 77kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human CAST

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

BS-17; PLACK;

Catalog #: A7634

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 831

Isotype: IgG

Swiss Prot: P20810

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is an endogenous calpain (calcium-dependent cysteine protease) inhibitor. It consists of an N-terminal domain L and four repetitive calpain-inhibition domains (domains 1-4), and it is involved in the proteolysis of amyloid precursor protein. The calpain/calpastatin system is involved in numerous membrane fusion events, such as neural vesicle exocytosis and platelet and red-cell aggregation. The encoded protein is also thought to affect the expression levels of genes encoding structural or regulatory proteins. Alternatively spliced transcript variants encoding different isoforms have been described.

To place an order, please [Click HERE](#).