

SLC2A4

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 55kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human SLC2A4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GLUT4;

Catalog #: A7637

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6517

Isotype: IgG

Swiss Prot: P14672

Purity: Affinity purification

For research use only.

Background:

This gene is a member of the solute carrier family 2 (facilitated glucose transporter) family and encodes a protein that functions as an insulin-regulated facilitative glucose transporter. In the absence of insulin, this integral membrane protein is sequestered within the cells of muscle and adipose tissue. Within minutes of insulin stimulation, the protein moves to the cell surface and begins to transport glucose across the cell membrane. Mutations in this gene have been associated with noninsulin-dependent diabetes mellitus (NIDDM).

To place an order, please [Click HERE](#).