

ANXA8L2

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 37kDa

Observed MW: Refer to figures

Immunogen:

Recombinant Protein of human ANXA8L2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ANXA8; ANXA8L2; VAC-beta; bA145E20.2;

Catalog #: A7641

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 728113

Isotype: IgG

Swiss Prot: Q5VT79

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the annexin family of evolutionarily conserved Ca²⁺ and phospholipid binding proteins. The encoded protein may function as an anticoagulant that indirectly inhibits the thromboplastin-specific complex. Overexpression of this gene has been associated with acute myelocytic leukemia. A highly similar duplicated copy of this gene is found in close proximity on the long arm of chromosome 10.

To place an order, please [Click HERE](#).