

CD7

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 25kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CD7

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GP40; TP41; Tp40; LEU-9;

Catalog #: A7650

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 924

Isotype: IgG

Swiss Prot: P09564

Purity: Affinity purification

For research use only.

Background:

This gene encodes a transmembrane protein which is a member of the immunoglobulin superfamily. This protein is found on thymocytes and mature T cells. It plays an essential role in T-cell interactions and also in T-cell/B-cell interaction during early lymphoid development.

To place an order, please [Click HERE](#).