

CHI3L2

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 44kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CHI3L2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CHIL2; YKL39; YKL-39;

Catalog #: A7654

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1117

Isotype: IgG

Swiss Prot: Q15782

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is similar to bacterial chitinases but lacks chitinase activity. The encoded protein is secreted and is involved in cartilage biogenesis. Several transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).