

CSF3R

Reactivity: Human Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 92kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human CSF3R

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CD114; GCSFR;

Catalog #: A7661

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1441

Isotype: IgG

Swiss Prot: Q99062

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is the receptor for colony stimulating factor 3, a cytokine that controls the production, differentiation, and function of granulocytes. The encoded protein, which is a member of the family of cytokine receptors, may also function in some cell surface adhesion or recognition processes. Alternatively spliced transcript variants have been described. Mutations in this gene are a cause of Kostmann syndrome, also known as severe congenital neutropenia.

To place an order, please [Click HERE](#).