

ARID3A

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 63kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ARID3A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

1 mg/ml

Synonym:

DRIL1; DRIL3; BRIGHT; E2FBP1;

Catalog #: A7668

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1820

Isotype: IgG

Swiss Prot: Q99856

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the ARID (AT-rich interaction domain) family of DNA binding proteins. It was found by homology to the Drosophila dead ringer gene, which is important for normal embryogenesis. Other ARID family members have roles in embryonic patterning, cell lineage gene regulation, cell cycle control, transcriptional regulation, and possibly in chromatin structure modification.

To place an order, please [Click HERE](#).