

GNAI2

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 40kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human GNAI2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GIP; GNAI2B; H_LUCA15.1; H_LUCA16.1;

Catalog #: A7676

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2771

Isotype: IgG

Swiss Prot: P04899

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is an alpha subunit of guanine nucleotide binding proteins (G proteins). The encoded protein contains the guanine nucleotide binding site and is involved in the hormonal regulation of adenylate cyclase. Several transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).