

HOXD10

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 38kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human HOXD10

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HOX4; HOX4D; HOX4E; Hox-4.4;

Catalog #: A7686

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3236

Isotype: IgG

Swiss Prot: P28358

Purity: Affinity purification

For research use only.

Background:

This gene is a member of the Abd-B homeobox family and encodes a protein with a homeobox DNA-binding domain. It is included in a cluster of homeobox D genes located on chromosome 2. The encoded nuclear protein functions as a sequence-specific transcription factor that is expressed in the developing limb buds and is involved in differentiation and limb development. Mutations in this gene have been associated with Wilm's tumor and congenital vertical talus (also known as "rocker-bottom foot" deformity or congenital convex pes valgus) and/or a foot deformity resembling that seen in Charcot-Marie-Tooth disease.

To place an order, please [Click HERE](#).