

INPP4A

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 110kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human INPP4A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

INPP4; TVAS1;

Catalog #: A7691

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3631

Isotype: IgG

Swiss Prot: Q96PE3

Purity: Affinity purification

For research use only.

Background:

This gene encodes an Mg⁺⁺ independent enzyme that hydrolyzes the 4-position phosphate from the inositol ring of phosphatidylinositol 3,4-bisphosphate, inositol 1,3,4-trisphosphate, and inositol 3,4-bisphosphate. Multiple transcript variants encoding distinct isoforms have been described.

To place an order, please [Click HERE](#).