

ITGA9

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 114kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ITGA9

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

q

Synonym:

RLC; ITGA4L; ALPHA-RLC;

Catalog #: A7693

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3680

Isotype: IgG

Swiss Prot: Q13797

Purity: Affinity purification

For research use only.

Background:

This gene encodes an alpha integrin. Integrins are heterodimeric integral membrane glycoproteins composed of an alpha chain and a beta chain that mediate cell-cell and cell-matrix adhesion. The protein encoded by this gene, when bound to the beta 1 chain, forms an integrin that is a receptor for VCAM1, cytotoxin and osteopontin. Expression of this gene has been found to be upregulated in small cell lung cancers.

To place an order, please [Click HERE](#).