

MATN3

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 53kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MATN3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HOA; OS2; EDM5; DIPOA; OADIP;

Catalog #: A7700

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4148

Isotype: IgG

Swiss Prot: O15232

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of von Willebrand factor A domain containing protein family. This family of proteins is thought to be involved in the formation of filamentous networks in the extracellular matrices of various tissues. This protein contains two von Willebrand factor A domains; it is present in the cartilage extracellular matrix and has a role in the development and homeostasis of cartilage and bone. Mutations in this gene result in multiple epiphyseal dysplasia.

To place an order, please [Click HERE](#).