

YBX1

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:100 - 1:200

Calculated MW:36kDa

Observed MW:Refer to figures

Immunogen:

A synthetic peptide of human YBX1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

jk

Synonym:

YB1; BP-8; CSDB; DBPB; YB-1; CSDA2; NSEP1; NSEP-1; MDR-NF1;

To place an order, please [Click HERE](#).

Catalog #:A7704

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:4904

Isotype:IgG

Swiss Prot:P67809

Purity:Affinity purification

For research use only.