

PPIB

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 24kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human PPIB

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

OI9; CYPB; SCYLP; CYP-S1; HEL-S-39;

Catalog #: A7713

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5479

Isotype: IgG

Swiss Prot: P23284

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a cyclosporine-binding protein and is mainly located within the endoplasmic reticulum. It is associated with the secretory pathway and released in biological fluids. This protein can bind to cells derived from T- and B-lymphocytes, and may regulate cyclosporine A-mediated immunosuppression. Variants have been identified in this protein that give rise to recessive forms of osteogenesis imperfecta.

To place an order, please [Click HERE](#).