

PRKACG

Reactivity:Human

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW:40kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human PRKACG

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

KAPG; PKACg; BDPLT19;

Catalog #:A7715

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:5568

Isotype:IgG

Swiss Prot:P22612

Purity:Affinity purification

For research use only.

Background:

Cyclic AMP-dependent protein kinase (PKA) consists of two catalytic subunits and a regulatory subunit dimer. This gene encodes the gamma form of its catalytic subunit. The gene is intronless and is thought to be a retrotransposon derived from the gene for the alpha form of the PKA catalytic subunit.

To place an order, please [Click HERE](#).