

AKR1B10

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 36kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human AKR1B10

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HIS; HSI; ARL1; ARL-1; ALDRLn; AKR1B11; AKR1B12;

Catalog #: A7823

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 57016

Isotype: IgG

Swiss Prot: O60218

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the aldo/keto reductase superfamily, which consists of more than 40 known enzymes and proteins. This member can efficiently reduce aliphatic and aromatic aldehydes, and it is less active on hexoses. It is highly expressed in adrenal gland, small intestine, and colon, and may play an important role in liver carcinogenesis.

To place an order, please [Click HERE](#).