

VANGL2

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 60kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human VANGL2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

LPP1; LTAP; STB1; STBM; STBM1;

Catalog #: A7825

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 57216

Isotype: IgG

Swiss Prot: Q9ULK5

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a membrane protein involved in the regulation of planar cell polarity, especially in the stereociliary bundles of the cochlea. The encoded protein transmits directional signals to individual cells or groups of cells in epithelial sheets. This protein is also involved in the development of the neural plate.

To place an order, please [Click HERE](#).