

RAB17

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 23kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RAB17

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A7831

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 64284

Isotype: IgG

Swiss Prot: Q9H0T7

Purity: Affinity purification

For research use only.

Background:

The Rab subfamily of small GTPases plays an important role in the regulation of membrane trafficking. RAB17 is an epithelial cell-specific GTPase (Lutcke et al., 1993)

To place an order, please [Click HERE](#).