

TBL1XR1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 56kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TBL1XR1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

C21; DC42; IRA1; TBLR1;

Catalog #: A7834

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 79718

Isotype: IgG

Swiss Prot: Q9BZK7

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene has sequence similarity with members of the WD40 repeat-containing protein family. The WD40 group is a large family of proteins, which appear to have a regulatory function. It is believed that the WD40 repeats mediate protein-protein interactions and members of the family are involved in signal transduction, RNA processing, gene regulation, vesicular trafficking, cytoskeletal assembly and may play a role in the control of cytotypic differentiation.

To place an order, please [Click HERE](#).