

DHDDS

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 39kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human DHDDS

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DS; CIT; CPT; HDS; RP59;

Catalog #: A7835

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 79947

Isotype: IgG

Swiss Prot: Q86SQ9

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene catalyzes cis-prenyl chain elongation to produce the polyprenyl backbone of dolichol, a glycosyl carrier lipid required for the biosynthesis of several classes of glycoproteins. Mutations in this gene are associated with retinitis pigmentosa type 59. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.

To place an order, please [Click HERE](#).