

FERMT3

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 76kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human FERMT3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

URP2; KIND3; MIG-2; MIG2B; URP2SF; UNC112C;

Catalog #: A7839

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 83706

Isotype: IgG

Swiss Prot: Q86UX7

Purity: Affinity purification

For research use only.

Background:

Kindlins are a small family of proteins that mediate protein-protein interactions involved in integrin activation and thereby have a role in cell adhesion, migration, differentiation, and proliferation. The protein encoded by this gene has a key role in the regulation of hemostasis and thrombosis. This protein may also help maintain the membrane skeleton of erythrocytes. Mutations in this gene cause the autosomal recessive leukocyte adhesion deficiency syndrome-III (LAD-III). Alternative splicing results in multiple transcript variants encoding distinct isoforms.

To place an order, please [Click HERE](#).