

PLCD4

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 88kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human PLCD4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A7841

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 84812

Isotype: IgG

Swiss Prot: Q9BRC7

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the delta class of phospholipase C enzymes. Phospholipase C enzymes play a critical role in many cellular processes by hydrolyzing phosphatidylinositol 4,5-bisphosphate into two intracellular second messengers, inositol 1,4,5-trisphosphate and diacylglycerol. Expression of this gene may be a marker for cancer.

To place an order, please [Click HERE](#).