

TCEAL8

Reactivity:Human

Tested applications:WB IF

Recommended Dilution:WB 1:500 - 1:2000 IF 1:50 - 1:100

Calculated MW:14kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human TCEAL8

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

WEX3;

Catalog #:A7845

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:90843

Isotype:IgG

Swiss Prot:Q8IYN2

Purity:Affinity purification

For research use only.

Background:

This gene encodes a member of the transcription elongation factor A (SII)-like (TCEAL) gene family. Members of this family contain TFA domains and may function as nuclear phosphoproteins that modulate transcription in a promoter context-dependent manner. Multiple family members are located on the X chromosome. Alternative splicing results in multiple transcript variants encoding a single isoform.

To place an order, please [Click HERE](#).