

## HAUS8

---

**Reactivity:** Human

**Tested applications:** WB IHC

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:** 45kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human HAUS8

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

DGT4; HICE1; NY-SAR-48;

**Catalog #:** A7847

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 93323

**Isotype:** IgG

**Swiss Prot:** Q9BT25

**Purity:** Affinity purification

For research use only.

**Background:**

HAUS8 is 1 of 8 subunits of the 390-kD human augmin complex, or HAUS complex. The augmin complex was first identified in *Drosophila*, and its name comes from the Latin verb 'augmentare,' meaning 'to increase.' The augmin complex is a microtubule-binding complex involved in microtubule generation within the mitotic spindle and is vital to mitotic spindle assembly (Goshima et al., 2008)

*To place an order, please [Click HERE](#).*