

LRG1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 38kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human LRG1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

LRG; HMFT1766;

Catalog #: A7850

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 116844

Isotype: IgG

Swiss Prot: P02750

Purity: Affinity purification

For research use only.

Background:

The leucine-rich repeat (LRR) family of proteins, including LRG1, have been shown to be involved in protein-protein interaction, signal transduction, and cell adhesion and development. LRG1 is expressed during granulocyte differentiation (O'Donnell et al., 2002)

To place an order, please [Click HERE](#).