

GLYCTK

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 55kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human GLYCTK

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HBEBP2; HBEBP4; HBeAgBP4A;

Catalog #: A7852

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 132158

Isotype: IgG

Swiss Prot: Q8IVS8

Purity: Affinity purification

For research use only.

Background:

This locus encodes a member of the glycerate kinase type-2 family. The encoded enzyme catalyzes the phosphorylation of (R)-glycerate and may be involved in serine degradation and fructose metabolism. Decreased activity of the encoded enzyme may be associated with the disease D-glyceric aciduria. Alternatively spliced transcript variants have been described.

To place an order, please [Click HERE](#).