

NDUFS7

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 24kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human NDUFS7

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PSST; CI-20; MY017; CI-20KD;

Catalog #: A7862

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 374291

Isotype: IgG

Swiss Prot: O75251

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein that is a subunit of one of the complexes that forms the mitochondrial respiratory chain. This protein is one of over 40 subunits found in complex I, the nicotinamide adenine dinucleotide (NADH):ubiquinone oxidoreductase. This complex functions in the transfer of electrons from NADH to the respiratory chain, and ubiquinone is believed to be the immediate electron acceptor for the enzyme. Mutations in this gene cause Leigh syndrome due to mitochondrial complex I deficiency, a severe neurological disorder that results in bilaterally symmetrical necrotic lesions in subcortical brain regions.

To place an order, please [Click HERE](#).