

AMPD1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 90kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human AMPD1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MAD; MADA; MMDD;

Catalog #: A7876

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 270

Isotype: IgG

Swiss Prot: P23109

Purity: Affinity purification

For research use only.

Background:

Adenosine monophosphate deaminase 1 catalyzes the deamination of AMP to IMP in skeletal muscle and plays an important role in the purine nucleotide cycle. Two other genes have been identified, AMPD2 and AMPD3, for the liver- and erythrocyte-specific isoforms, respectively.

Deficiency of the muscle-specific enzyme is apparently a common cause of exercise-induced myopathy and probably the most common cause of metabolic myopathy in the human.

Alternatively spliced transcript variants encoding different isoforms have been identified in this gene.

To place an order, please [Click HERE](#).