

CHD1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:100 - 1:200

Calculated MW: 197kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CHD1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Catalog #: A7883

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1105

Isotype: IgG

Swiss Prot: O14646

Purity: Affinity purification

For research use only.

Background:

The CHD family of proteins is characterized by the presence of chromo (chromatin organization modifier) domains and SNF2-related helicase/ATPase domains. CHD genes alter gene expression possibly by modification of chromatin structure thus altering access of the transcriptional apparatus to its chromosomal DNA template.

To place an order, please [Click HERE](#).