

CLTC

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 192kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human CLTC

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

Hc; CHC; CHC17; CLH-17; CLTCL2;

Catalog #: A7886

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1213

Isotype: IgG

Swiss Prot: Q00610

Purity: Affinity purification

For research use only.

Background:

Clathrin is a major protein component of the cytoplasmic face of intracellular organelles, called coated vesicles and coated pits. These specialized organelles are involved in the intracellular trafficking of receptors and endocytosis of a variety of macromolecules. The basic subunit of the clathrin coat is composed of three heavy chains and three light chains.

To place an order, please [Click HERE](#).