

NCAM1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:100 - 1:200

Calculated MW: 95kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human NCAM1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

i

Synonym:

CD56; NCAM; MSK39;

Catalog #: A7913

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4684

Isotype: IgG

Swiss Prot: P13591

Purity: Affinity purification

For research use only.

Background:

This gene encodes a cell adhesion protein which is a member of the immunoglobulin superfamily. The encoded protein is involved in cell-to-cell interactions as well as cell-matrix interactions during development and differentiation. The encoded protein has been shown to be involved in development of the nervous system, and for cells involved in the expansion of T cells and dendritic cells which play an important role in immune surveillance. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).