

SOCS5

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 61kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SOCS5

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CIS6; CISH6; Cish5; SOCS-5;

Catalog #: A7952

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9655

Isotype: IgG

Swiss Prot: O75159

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene contains a SH2 domain and a SOCS BOX domain. The protein thus belongs to the suppressor of cytokine signaling (SOCS) family, also known as STAT-induced STAT inhibitor (SSI) protein family. SOCS family members are known to be cytokine-inducible negative regulators of cytokine signaling. The specific function of this protein has not yet been determined. Two alternatively spliced transcript variants encoding an identical protein have been reported.

To place an order, please [Click HERE](#).