

MTX2

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 30kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MTX2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

h

Catalog #: A7958

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10651

Isotype: IgG

Swiss Prot: O75431

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is highly similar to the metaxin 2 protein from mouse, which has been shown to interact with the mitochondrial membrane protein metaxin 1. Because of this similarity, it is thought that the encoded protein is peripherally associated with the cytosolic face of the outer mitochondrial membrane, and that it is involved in the import of proteins into the mitochondrion. Alternative splicing results in multiple transcript variants. A related pseudogene has been identified on chromosome 7.

To place an order, please [Click HERE](#).