

ASCC3

Reactivity: Human Mouse

Tested applications: WB IF

Recommended Dilution: WB 1:500 - 1:2000 IF 1:50 - 1:100

Calculated MW: 251kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ASCC3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

n

Synonym:

RNAH; HELIC1; ASC1p200;

Catalog #: A7960

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10973

Isotype: IgG

Swiss Prot: Q8N3C0

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein that belongs to a family of helicases that are involved in the ATP-dependent unwinding of nucleic acid duplexes. The encoded protein is the largest subunit of the activating signal cointegrator 1 complex that is involved in DNA repair and resistance to alkylation damage. Alternate splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).