

PRDM7

Reactivity: Human Mouse Rat

Tested applications: WB IF

Recommended Dilution: WB 1:500 - 1:2000 IF 1:50 - 1:100

Calculated MW: 56kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human PRDM7

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PFM4; ZNF910;

Catalog #: A7961

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 11105

Isotype: IgG

Swiss Prot: Q9NQW5

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of a family of proteins that may have roles in transcription and other nuclear processes. The encoded protein contains a KRAB (Kruppel-associated box) domain -A box and a SET (Su(var)3-9, Enhancer-of-zeste, Trithorax) domain and may function as a histone methyltransferase.

To place an order, please [Click HERE](#).