

B9D1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 23kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human B9D1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

B9; MKS9; EPPB9; MKSR1;

Catalog #: A7973

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 27077

Isotype: IgG

Swiss Prot: Q9UPM9

Purity: Affinity purification

For research use only.

Background:

This gene encodes a B9 domain-containing protein, one of several that are involved in ciliogenesis. Alterations in expression of this gene have been found in a family with Meckel syndrome. Meckel syndrome has been associated with at least six different genes. This gene is located within the Smith-Magenis syndrome region on chromosome 17. Three alternatively spliced transcript variants that encode different proteins have been described for this gene.

To place an order, please [Click HERE](#).