

## USP25

---

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:** 122kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human USP25

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

USP21;

**Catalog #:** A7975

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 29761

**Isotype:** IgG

**Swiss Prot:** Q9UHP3

**Purity:** Affinity purification

For research use only.

**Background:**

Ubiquitin (MIM 191339) is a highly conserved 76-amino acid protein involved in regulation of intracellular protein breakdown, cell cycle regulation, and stress response. Ubiquitin is released from degraded proteins by disassembly of the polyubiquitin chains, which is mediated by ubiquitin-specific proteases (USPs), such as USP25 (Valero et al., 1999)

*To place an order, please [Click HERE](#).*