

ZDHHC9

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 41kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ZDHHC9

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CGI89; DHHC9; MMSA1; MRXSZ; ZNF379; ZNF380; CXorf11; ZDHHC10;

Catalog #: A7977

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 51114

Isotype: IgG

Swiss Prot: Q9Y397

Purity: Affinity purification

For research use only.

Background:

This gene encodes an integral membrane protein that is a member of the zinc finger DHHC domain-containing protein family. The encoded protein forms a complex with golgin subfamily A member 7 and functions as a palmitoyltransferase. This protein specifically palmitoylates HRAS and NRAS. Mutations in this gene are associated with X-linked mental retardation. Alternate splicing results in multiple transcript variants that encode the same protein.

To place an order, please [Click HERE](#).