

TERF2IP

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 44kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TERF2IP

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RAP1; DRIP5;

Catalog #: A7981

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 54386

Isotype: IgG

Swiss Prot: Q9NYB0

Purity: Affinity purification

For research use only.

Background:

The gene encodes a protein that is part of a complex involved in telomere length regulation. Pseudogenes are present on chromosomes 5 and 22.

To place an order, please [Click HERE](#).