

EGLN3

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 27kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human EGLN3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

bf

Synonym:

PHD3; HIFPH3; HIFP4H3;

To place an order, please [Click HERE](#).

Catalog #: A8001

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 112399

Isotype: IgG

Swiss Prot: Q9H6Z9

Purity: Affinity purification

For research use only.