

SOCS4

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:51kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human SOCS4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SOCS7;

Catalog #:A8003

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:122809

Isotype:IgG

Swiss Prot:Q8WXH5

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene contains a SH2 domain and a SOCS BOX domain. The protein thus belongs to the suppressor of cytokine signaling (SOCS), also known as STAT-induced STAT inhibitor (SSI), protein family. SOCS family members are known to be cytokine-inducible negative regulators of cytokine signaling. Two alternatively spliced transcript variants encoding the same protein have been found for this gene.

To place an order, please [Click HERE](#).