

GRID1

Reactivity: Human Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 112kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human GRID1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GluD1;

Catalog #: A8011

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2894

Isotype: IgG

Swiss Prot: Q9ULK0

Purity: Affinity purification

For research use only.

Background:

This gene encodes a subunit of glutamate receptor channels. These channels mediate most of the fast excitatory synaptic transmission in the central nervous system and play key roles in synaptic plasticity.

To place an order, please [Click HERE](#).