

PRICKLE2

Reactivity: Human

Tested applications: WB IF

Recommended Dilution: WB 1:500 - 1:2000 IF 1:50 - 1:100

Calculated MW: 96kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human PRICKLE2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

EPM5;

Catalog #: A8023

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 166336

Isotype: IgG

Swiss Prot: Q7Z3G6

Purity: Affinity purification

For research use only.

Background:

This gene encodes a homolog of Drosophila prickle. The exact function of this gene is not known, however, studies in mice suggest that it may be involved in seizure prevention. Mutations in this gene are associated with progressive myoclonic epilepsy type 5.

To place an order, please [Click HERE](#).