

## DIS3

---

**Reactivity:**Human Rat

**Tested applications:**WB

**Recommended Dilution:**WB 1:1000 - 1:2000

**Calculated MW:**109kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human DIS3

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

RRP44; dis3p; EXOSC11; KIAA1008; 2810028N01Rik;

*To place an order, please [Click HERE](#).*

**Catalog #:**A8027

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**22894

**Isotype:**IgG

**Swiss Prot:**Q9Y2L1

**Purity:**Affinity purification

For research use only.