

RALYL

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:200 - 1:2000 IHC 1:20 - 1:200

Calculated MW: 32kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human RALYL

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HNRPCL3;

To place an order, please [Click HERE](#).

Catalog #: A8032

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 138046

Isotype: IgG

Swiss Prot: Q86SE5

Purity: Affinity purification

For research use only.