

ENPP3

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:1000

Calculated MW: 100kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human ENPP3

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

B10; NPP3; PDNP3; CD203c; PD-IBETA;

Catalog #: A8042

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5169

Isotype: IgG

Swiss Prot: O14638

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene belongs to a series of ectoenzymes that are involved in hydrolysis of extracellular nucleotides. These ectoenzymes possess ATPase and ATP pyrophosphatase activities and are type II transmembrane proteins. Expression of the related rat mRNA has been found in a subset of immature glial cells and in the alimentary tract. The corresponding rat protein has been detected in the pancreas, small intestine, colon, and liver. The human mRNA is expressed in glioma cells, prostate, and uterus. Expression of the human protein has been detected in uterus, basophils, and mast cells.

To place an order, please [Click HERE](#).