

LPO

Reactivity:Human Mouse

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:80kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human LPO

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SPO;

Catalog #:A8053

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:4025

Isotype:IgG

Swiss Prot:P22079

Purity:Affinity purification

For research use only.

Background:

This gene encodes an oxidoreductase secreted from salivary, mammary, and other mucosal glands that functions as a natural antibacterial agent. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).