

MED24

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 110kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MED24

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CRSP4; ARC100; THRAP4; CRSP100; DRIP100; TRAP100;

Catalog #: A8059

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9862

Isotype: IgG

Swiss Prot: O75448

Purity: Affinity purification

For research use only.

Background:

This gene encodes a component of the mediator complex (also known as TRAP, SMCC, DRIP, or ARC), a transcriptional coactivator complex thought to be required for the expression of almost all genes. The mediator complex is recruited by transcriptional activators or nuclear receptors to induce gene expression, possibly by interacting with RNA polymerase II and promoting the formation of a transcriptional pre-initiation complex. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).